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Instructions for use

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Positioning Comfortable Plus Duo
Velcro Stabilo System

Typ

Comfortable Plus Duo Velcro
inna nazwa: Stairclimber Duo Velcro
Comfortable Plus Duo Hilo Velcro

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Thank you for choosing our product. The NOEL Żebrowski sp. j. has designed and manufactured Comfortable Plus Duo Velcro Stabulo System to the highest safety and functional use standards. We are confident that our product will meet your highest requirements and expectations. These instructions for use are an integral part of vacuum positioning system stabilizing torso, limbs and head Comfortable Plus Duo Velcro Stabulo System. Please read the instructions thoroughly before using the product.

Vacuum positioning system stabilizing torso, limbs and head Comfortable Plus Duo Velcro Stabulo System

Vacuum Positioning Comfortable Plus Duo Velcro Stabulo System is a product designed for children and adults with cerebral palsy and myelomeningocele, compensating for and alleviating the disability by stabilising torso, head and limbs in sitting position. Positioning Comfortable Plus Duo Velcro Stabulo System can be used for compensating torso distortions and stabilizing the flaccid patient.

INDICATIONS

- cerebral palsy
- myelomeningocele
- torso distortions
- flaccidity

Comfortable Plus Duo Velcro Stabulo System cushions are covered a special outer material that allows attaching velcro strips. The Velcro strap is part of the Comfortable Plus Duo Stabulo System. Comfortable Plus Duo Velcro Hilo Stabulo System has an additional abduction wedge that reduces contractures and ensures proper placement of the legs in a sitting position.

PRODUCT SAFETY WARNINGS

- positioning Comfortable Plus Duo Velcro Stabulo System must be securely fixed to seats or backrests
- use of airy cover is advisable if increased ventilation is required around the body
- avoid contact with sharp objects that may damage the cushion
- do not clean using solvent-based cleaners, e.g. toluene or acetone
- do not pull on valve
- do not overinflate cushion
- use only special pump included with cushion for inflation

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- do not iron
- in case of tear, micro-granules may escape which may be dangerous for small children
- only persons who received instructions or training may use the device

CONTRAINDICATIONS

Permitted user weight 80 kg

ATTENTION

- the thickness of the cushion can naturally reduce by 5% in the initial period of use
- in case of considerable changes of temperature use pump again to obtain proper hardness (higher temperature softens the cushion, lower temperature hardens the cushion)

INSTRUCTIONS FOR USE

Hand pump

1. Take out tube in pump handle.
2. Screw tube into cushion valve and wait until cushion reaches the required softness.
3. Shape cushion (either manually to adjust shape or by sitting or leaning on it).
4. Screw pump into tube and deflate air until desired hardness has been reached.
5. Remove tube from pump and screw protecting cap on valve.
6. To change shape - screw tube into valve and wait until cushion softens, then repeat steps from 2 to 5.

Electric pump

1. Take out tube from box with electric pump.
2. Screw off valve cap, screw tube (tapped part) into cushion valve, and push second tube into SOFT socket in electric pump. Turn on pump and wait until cushion will be soft enough, then turn pump off.
3. Start shaping cushion (either manually to adjust shape or by sitting or leaning on it).
4. Change position of tube from SOFT to HARD and turn pump on removing air until the desired hardness has been reached.
5. Screw off tube, put cap back on and tighten.
6. To change the shape repeat steps from 2 to 5.

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MAINTENANCE

- regularly check condition of product – pay attention to any cracks and rips
- screw on valve cap before washing
- hand wash cushion under running water using a standard detergent e.g. a washing up liquid
- suitable for washing in washing machine at 40°C cycle. Fabric softener is sufficient, does not require use of washing powder
- wipe dry with paper towel
- standard disinfectants used in hospitals can be used
- do not use solvent-based cleaners such as toluene or acetone
- for Velcro products, Hoover the soft surface regularly to prevent loss of adherence

WARRANTY

- the warranty is provided for a period of 24 months
- the warranty does not cover mechanical damage to Positioning Comfortable Plus Velcro Duo Stabilo System
- if in doubt, please contact seller or manufacturer directly
- report any incidents to seller or manufacturer directly

DECLARATION OF CONFORMITY

Vacuum positioning system stabilizing torso, limbs and head Comfortable Plus Duo Velcro Stabilo System medical device complies with the basic requirements of the Regulation (EU) 2017/745 of the European Parliament and of the Council of 5 April 2017 on medical devices.

NOEL Żebrowski sp. j.
ul. Christo Botewa 4A, 30-798 Kraków
+48 12 341 53 16
www.stabilo.krakow.pl
stabilo@stabilo.krakow.pl

Find Stabilo System Cushions on:

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